

Product datasheet for **GA211230**

Mouse *Pik3ap1* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3ap1 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase adaptor protein 1
Symbol:	Pik3ap1
UniProt ID:	Q9EQ32
Synonyms:	1810044J04Rik; BCAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)