

## Product datasheet for **GA211248**

### Mouse Cnnm1 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Cnnm1 CRISPRa kit - CRISPR gene activation of mouse cyclin M1                                                                                                                                                                                                                                              |
| Symbol:       | Cnnm1                                                                                                                                                                                                                                                                                                      |
| Components:   | GA211248G1, Cnnm1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211248G2, Cnnm1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211248G3, Cnnm1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q0GA42</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | Acdp1; AI846314; CLP-1                                                                                                                                                                                                                                                                                     |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)