

Product datasheet for **GA211250**

Mouse Pde4dip activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pde4dip CRISPRa kit - CRISPR gene activation of mouse phosphodiesterase 4D interacting protein (myomegalin)
Symbol:	Pde4dip
UniProt ID:	Q80YT7
Synonyms:	4732458A06Rik; 9430063L05Rik; C87016; D3Bwg1078e; D130016K21Rik; mKIAA0454; Usmg4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)