

Product datasheet for **GA211282**

Mouse Pi4k2a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pi4k2a CRISPRa kit - CRISPR gene activation of mouse phosphatidylinositol 4-kinase type 2 alpha
Symbol:	Pi4k2a
UniProt ID:	<u>Q2TBE6</u>
Synonyms:	Pi4k2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)