

Product datasheet for **GA211300**

Mouse Tac4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tac4 CRISPRa kit - CRISPR gene activation of mouse tachykinin 4
Symbol:	Tac4
Components:	GA211300G1, Tac4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211300G2, Tac4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211300G3, Tac4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99N14</u>
Synonyms:	AW489379; PPT-C
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)