

Product datasheet for GA211313

Mouse Lmod1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Lmod1 CRISPRa kit - CRISPR gene activation of mouse leiomodin 1 (smooth muscle)
Symbol:	Lmod1
UniProt ID:	Q8BVA4
Synonyms:	9530015K06Rik; SM-Lmod
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)