

Product datasheet for **GA211343**

Mouse Mpv17l activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Mpv17l CRISPRa kit - CRISPR gene activation of mouse Mpv17 transgene, kidney disease mutant-like
Symbol:	Mpv17l
Components:	GA211343G1, Mpv17l gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211343G2, Mpv17l gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211343G3, Mpv17l gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99MS3</u>
Synonyms:	M-LP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)