


Product datasheet for **GA211368**

Mouse Pcdhb3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcdhb3 CRISPRa kit - CRISPR gene activation of mouse protocadherin beta 3
Symbol:	Pcdhb3
Components:	GA211368G1, Pcdhb3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211368G2, Pcdhb3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211368G3, Pcdhb3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	PcdhbC
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)