

Product datasheet for **GA211374**

Mouse Pcdhb10 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcdhb10 CRISPRa kit - CRISPR gene activation of mouse protocadherin beta 10
Symbol:	Pcdhb10
Components:	GA211374G1, Pcdhb10 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211374G2, Pcdhb10 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211374G3, Pcdhb10 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Pcdhb5D; PcdhbJ
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)