

Product datasheet for **GA211378**

Mouse Pcdhb14 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcdhb14 CRISPRa kit - CRISPR gene activation of mouse protocadherin beta 14
Symbol:	Pcdhb14
Components:	GA211378G1, Pcdhb14 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211378G2, Pcdhb14 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211378G3, Pcdhb14 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q6PB90</u>
Synonyms:	2210006M07Rik; Pcdhb17; PcdhbN
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)