

Product datasheet for **GA211384**

Mouse Pcdhb20 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pcdhb20 CRISPRa kit - CRISPR gene activation of mouse protocadherin beta 20
Symbol:	Pcdhb20
Components:	GA211384G1, Pcdhb20 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211384G2, Pcdhb20 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211384G3, Pcdhb20 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
Synonyms:	Pcdhb14; PcdhbT
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)