

Product datasheet for **GA211410**

Mouse Trim34a activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Trim34a CRISPRa kit - CRISPR gene activation of mouse tripartite motif-containing 34A
Symbol:	Trim34a
Components:	GA211410G1, Trim34a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211410G2, Trim34a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211410G3, Trim34a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q99PP6</u>
Synonyms:	Trim34; Trim34-1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)