

Product datasheet for **GA211435**

Mouse Dgcr8 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dgcr8 CRISPRa kit - CRISPR gene activation of mouse DGCR8, microprocessor complex subunit
Symbol:	Dgcr8
UniProt ID:	Q9EQM6
Synonyms:	D16H22S788E; D16H22S1742E; D16Wis2; Gy1; mir-1306; N41; Vo59c07
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)