


Product datasheet for **GA211462**

Mouse Tprn activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tprn CRISPRa kit - CRISPR gene activation of mouse taperin
Symbol:	Tprn
Components:	GA211462G1, Tprn gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211462G2, Tprn gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211462G3, Tprn gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2AI08</u>
Synonyms:	C87750; C430004E15Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)