

Product datasheet for **GA211467**

Mouse C77080 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	C77080 CRISPRa kit - CRISPR gene activation of mouse expressed sequence C77080
Symbol:	C77080
Components:	GA211467G1, C77080 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211467G2, C77080 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211467G3, C77080 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>A2A7S8</u>
Synonyms:	Kiaa1522; mKIAA1522
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)