

Product datasheet for **GA211497**

Mouse D630023F18Rik activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	D630023F18Rik CRISPRa kit - CRISPR gene activation of mouse RIKEN cDNA D630023F18 gene
Symbol:	D630023F18Rik
Components:	GA211497G1, D630023F18Rik gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211497G2, D630023F18Rik gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211497G3, D630023F18Rik gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8C3M9</u>
Synonyms:	AI314969
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)