


Product datasheet for **GA211506**

Mouse Sh3bp4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Sh3bp4 CRISPRa kit - CRISPR gene activation of mouse SH3-domain binding protein 4
Symbol:	Sh3bp4
Components:	GA211506G1, Sh3bp4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211506G2, Sh3bp4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211506G3, Sh3bp4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q921I6
Synonyms:	AI594717; AW227605; BOG25
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)