


Product datasheet for **GA211522**

Mouse *Kcnc2* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnc2</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage gated channel, Shab-related subfamily, member 2
Symbol:	<i>Kcnc2</i>
Components:	GA211522G1, <i>Kcnc2</i> gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211522G2, <i>Kcnc2</i> gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211522G3, <i>Kcnc2</i> gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A6H8H5
Synonyms:	9630047L19Rik; BB130875; Kv2.2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)