

Product datasheet for GA211544

Mouse Hrh3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Hrh3 CRISPRa kit - CRISPR gene activation of mouse histamine receptor H3
Symbol:	Hrh3
UniProt ID:	P58406
Synonyms:	AW049250; H3R; HH3R
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)