


Product datasheet for **GA211555**

Mouse Wdr47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Wdr47 CRISPRa kit - CRISPR gene activation of mouse WD repeat domain 47
Symbol:	Wdr47
Components:	GA211555G1, Wdr47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211555G2, Wdr47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211555G3, Wdr47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CGF6</u>
Synonyms:	1810073M12Rik; AA552998; mKIAA0893
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)