

Product datasheet for **GA211568**

Mouse *Kcnc4* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kcnc4</i> CRISPRa kit - CRISPR gene activation of mouse potassium voltage gated channel, Shaw-related subfamily, member 4
Symbol:	<i>Kcnc4</i>
UniProt ID:	Q8R1C0
Synonyms:	A1850292; Kcr2-4; KSHIIC; Kv3.4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

