

Product datasheet for **GA211572**

Mouse *Ifi44* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ifi44 CRISPRa kit - CRISPR gene activation of mouse interferon-induced protein 44
Symbol:	Ifi44
UniProt ID:	Q8BV66
Synonyms:	A430056A10Rik; AW261460; MTAP44; p44
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)