


## Product datasheet for **GA211587**

### Mouse Phactr4 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                                  |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                                     |
| Description:  | Phactr4 CRISPRa kit - CRISPR gene activation of mouse phosphatase and actin regulator 4                                                                                                                                                                                                                          |
| Symbol:       | Phactr4                                                                                                                                                                                                                                                                                                          |
| Components:   | GA211587G1, Phactr4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211587G2, Phactr4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211587G3, Phactr4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q501J7</a>                                                                                                                                                                                                                                                                                           |
| Synonyms:     | 3110001B12Rik; AI527228; AW495572; C330013F19Rik; mKIAA4120; N28169                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                          |

#### Product images:



[View online »](#)