

Product datasheet for **GA211609**

Mouse Upk3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Upk3b CRISPRa kit - CRISPR gene activation of mouse uroplakin 3B
Symbol:	Upk3b
Components:	GA211609G1, Upk3b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211609G2, Upk3b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211609G3, Upk3b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q80YF6</u>
Synonyms:	AI447729; AI957187; P35; PMS2L14; UpIIIb
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)