

Product datasheet for **GA211632**

Mouse Ttl3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ttl3 CRISPRa kit - CRISPR gene activation of mouse tubulin tyrosine ligase-like family, member 3
Symbol:	Ttl3
Components:	GA211632G1, Ttl3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211632G2, Ttl3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211632G3, Ttl3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	A4Q9E5
Synonyms:	4833441J24Rik; AI450050
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)