

Product datasheet for **GA211633**

Mouse Tmem168 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Tmem168 CRISPRa kit - CRISPR gene activation of mouse transmembrane protein 168
Symbol:	Tmem168
Components:	GA211633G1, Tmem168 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211633G2, Tmem168 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211633G3, Tmem168 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91VX9</u>
Synonyms:	5730526F17Rik; 8430437G11Rik; AI462344; AI504145
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)