

Product datasheet for **GA211637**

Mouse Parp11 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Parp11 CRISPRa kit - CRISPR gene activation of mouse poly (ADP-ribose) polymerase family, member 11
Symbol:	Parp11
Components:	GA211637G1, Parp11 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211637G2, Parp11 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211637G3, Parp11 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8CFF0</u>
Synonyms:	5330431N24Rik; AI851877; ARTD11; HIN1L; PARP-11
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)