


Product datasheet for **GA211643**

Mouse Prrt4 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Prrt4 CRISPRa kit - CRISPR gene activation of mouse proline-rich transmembrane protein 4
Symbol:	Prrt4
Components:	GA211643G1, Prrt4 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211643G2, Prrt4 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211643G3, Prrt4 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RU40</u>
Synonyms:	AW455481; D330027H18Rik
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)