

Product datasheet for **GA211677**

Mouse Sf3b3 activation kit by CRISPRa

Product data:

| | |
|---------------|-----------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits |
| Description: | Sf3b3 CRISPRa kit - CRISPR gene activation of mouse splicing factor 3b, subunit 3 |
| Symbol: | Sf3b3 |
| UniProt ID: | <u>Q921M3</u> |
| Synonyms: | 1810061H24Rik; 5730409A01Rik; AA409318; D8Ertd633e; mKIAA0017; RSE1; SAP130 |
| Format: | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

Product images:



[View online >](#)