


## Product datasheet for **GA211681**

### Mouse *Ces2a* activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Ces2a CRISPRa kit - CRISPR gene activation of mouse carboxylesterase 2A                                                                                                                                                                                                                                    |
| Symbol:       | Ces2a                                                                                                                                                                                                                                                                                                      |
| Components:   | GA211681G1, Ces2a gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211681G2, Ces2a gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211681G3, Ces2a gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <u><a href="#">Q8QZR3</a></u>                                                                                                                                                                                                                                                                              |
| Synonyms:     | 9130231C15Rik; AI266984; Ces6                                                                                                                                                                                                                                                                              |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)