


Product datasheet for **GA211697**

Mouse Cyp4v3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Cyp4v3 CRISPRa kit - CRISPR gene activation of mouse cytochrome P450, family 4, subfamily v, polypeptide 3
Symbol:	Cyp4v3
Components:	GA211697G1, Cyp4v3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211697G2, Cyp4v3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211697G3, Cyp4v3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q9DBW0
Synonyms:	AU043077; AW111961; Cyp4v2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)