

Product datasheet for **GA211706**

Mouse Pls1 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Pls1 CRISPRa kit - CRISPR gene activation of mouse plastin 1 (I-isoform)
Symbol:	Pls1
Components:	GA211706G1, Pls1 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211706G2, Pls1 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211706G3, Pls1 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q3V0K9</u>
Synonyms:	AI427122
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)