

Product datasheet for **GA211717**

Mouse **Lrrc49** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Lrrc49 CRISPRa kit - CRISPR gene activation of mouse leucine rich repeat containing 49
Symbol:	Lrrc49
Components:	GA211717G1, Lrrc49 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211717G2, Lrrc49 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211717G3, Lrrc49 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q91YK0</u>
Synonyms:	AI430817; AW050057; D430025H09Rik; p79; PGs4
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)