

## Product datasheet for GA211721

### Mouse Pls3 activation kit by CRISPRa

#### Product data:

|               |   |
|---------------|---|
| Product Type: | CRISPRa Kits  |
| Description:  | Pls3 CRISPRa kit - CRISPR gene activation of mouse platin 3 (T-isoform) |
| Symbol:       | Pls3  |
| UniProt ID:   | <a href="#">Q99K51</a>  |
| Synonyms:     | A1115446; AL024105  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug) |

#### Product images:



[View online >](#)