

## Product datasheet for **GA211738**

### Mouse Fig4 activation kit by CRISPRa

#### Product data:

|               |  |
|---------------|--|
| Product Type: | CRISPRa Kits   |
| Description:  | Fig4 CRISPRa kit - CRISPR gene activation of mouse FIG4 phosphoinositide 5-phosphatase |
| Symbol:       | Fig4   |
| UniProt ID:   | <u><a href="#">Q91WF7</a></u>  |
| Synonyms:     | A530089117Rik; AI326867; Sac3  |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                |

#### Product images:



[View online >](#)