

## Product datasheet for **GA211741**

### Mouse *Ttc41* activation kit by CRISPRa

#### Product data:

|               |                                                                                        |
|---------------|----------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                           |
| Description:  | Ttc41 CRISPRa kit - CRISPR gene activation of mouse tetratricopeptide repeat domain 41 |
| Symbol:       | Ttc41                                                                                  |
| UniProt ID:   | <a href="#">Q692V3</a>                                                                 |
| Synonyms:     | AI449705; AU041509; Gnn                                                                |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                |

#### Product images:



[View online >](#)