

Product datasheet for GA211747

Mouse Nt5dc3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Nt5dc3 CRISPRa kit – CRISPR gene activation of mouse 5'-nucleotidase domain containing 3
Symbol:	Nt5dc3
Components:	GA211747G1, Nt5dc3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211747G2, Nt5dc3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211747G3, Nt5dc3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UHB1
Synonyms:	AU040402; AW540062; C630002B14Rik; TUI2B1-TY
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

