

Product datasheet for GA211763

Mouse Slc25a41 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a41 CRISPRa kit – CRISPR gene activation of mouse solute carrier family 25, member 41
Symbol:	Slc25a41
Components:	GA211763G1, Slc25a41 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211763G2, Slc25a41 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211763G3, Slc25a41 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8BVN7</u>
Synonyms:	4933406J04Rik; AI507137
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:

