

## Product datasheet for **GA211818**

### Mouse Rexo2 activation kit by CRISPRa

#### Product data:

|               |                                                                                                                                                                                                                                                                                                            |
|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Type: | CRISPRa Kits                                                                                                                                                                                                                                                                                               |
| Description:  | Rexo2 CRISPRa kit - CRISPR gene activation of mouse RNA exonuclease 2                                                                                                                                                                                                                                      |
| Symbol:       | Rexo2                                                                                                                                                                                                                                                                                                      |
| Components:   | GA211818G1, Rexo2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211818G2, Rexo2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074)<br>GA211818G3, Rexo2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074)<br>1 CRISPRa-Enhancer vector, SKU GE100056<br>1 CRISPRa scramble vector, SKU GE100077 |
| UniProt ID:   | <a href="#">Q9D8S4</a>                                                                                                                                                                                                                                                                                     |
| Synonyms:     | 1810038D15Rik; AW107347; Rex2; Sfn; Smfn                                                                                                                                                                                                                                                                   |
| Format:       | 3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)                                                                                                                                                                                                                                    |

#### Product images:



[View online »](#)