


Product datasheet for **GA211823**

Mouse Ppp4r3b activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ppp4r3b CRISPRa kit - CRISPR gene activation of mouse protein phosphatase 4 regulatory subunit 3B
Symbol:	Ppp4r3b
Components:	GA211823G1, Ppp4r3b gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211823G2, Ppp4r3b gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211823G3, Ppp4r3b gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q922R5</u>
Synonyms:	AW011752; AW557776; mKIAA1387; Smek2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)