

Product datasheet for **GA211829**

Mouse **Pik3r6** activation kit by **CRISPRa**

Product data:

Product Type:	CRISPRa Kits
Description:	Pik3r6 CRISPRa kit - CRISPR gene activation of mouse phosphoinositide-3-kinase regulatory subunit 5
Symbol:	Pik3r6
Components:	GA211829G1, Pik3r6 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211829G2, Pik3r6 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211829G3, Pik3r6 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3U6Q4
Synonyms:	BB220380; p84; p84 PIKAP; p87(PIKAP); p87PIKAP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)