

Product datasheet for **GA211838**

Mouse Ptpn23 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Ptpn23 CRISPRa kit - CRISPR gene activation of mouse protein tyrosine phosphatase, non-receptor type 23
Symbol:	Ptpn23
UniProt ID:	Q6PB44
Synonyms:	A1462446; PTP-TD14
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)