

Product datasheet for **GA211845**

Mouse Slc25a47 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Slc25a47 CRISPRa kit - CRISPR gene activation of mouse solute carrier family 25, member 47
Symbol:	Slc25a47
Components:	GA211845G1, Slc25a47 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211845G2, Slc25a47 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211845G3, Slc25a47 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q6IS41
Synonyms:	AI132487; AI876593; HDMCP
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)