

Product datasheet for **GA211855**

Mouse Rnf44 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rnf44 CRISPRa kit - CRISPR gene activation of mouse ring finger protein 44
Symbol:	Rnf44
Components:	GA211855G1, Rnf44 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211855G2, Rnf44 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211855G3, Rnf44 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	Q3UHI8
Synonyms:	AI854545; mKIAA1100
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)