

Product datasheet for **GA211893**

Mouse *Kdelr3* activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	<i>Kdelr3</i> CRISPRa kit - CRISPR gene activation of mouse KDEL (Lys-Asp-Glu-Leu) endoplasmic reticulum protein retention receptor 3
Symbol:	<i>Kdelr3</i>
UniProt ID:	Q8R1L4
Synonyms:	A1173274
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)