

## Product datasheet for **GA211911**

### Mouse Topors activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Topors CRISPRa kit - CRISPR gene activation of mouse topoisomerase I binding, arginine/serine-rich
Symbol:	Topors
Components:	GA211911G1, Topors gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211911G2, Topors gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211911G3, Topors gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u><a href="#">Q80Z37</a></u>
Synonyms:	AW105885; LUN; p53BP3/LUN; TP53BPL
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:



[View online »](#)