

Product datasheet for **GA211923**

Mouse Rrn3 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Rrn3 CRISPRa kit - CRISPR gene activation of mouse RRN3 RNA polymerase I transcription factor homolog (yeast)
Symbol:	Rrn3
Components:	GA211923G1, Rrn3 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211923G2, Rrn3 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211923G3, Rrn3 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>B2RS91</u>
Synonyms:	AL023001; E130302O19Rik; R75565; TIF-1A; Tif1a
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)