

Product datasheet for GA211930

Mouse Stk38 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Stk38 CRISPRa kit - CRISPR gene activation of mouse serine/threonine kinase 38
Symbol:	Stk38
UniProt ID:	Q91VJ4
Synonyms:	5830476G13Rik; 9530097A09Rik; AA617404; Ndr1
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online >](#)