

## Product datasheet for **GA211933**

### Mouse Tocr activation kit by CRISPRa

#### Product data:

Product Type:	CRISPRa Kits
Description:	Tocr CRISPRa kit - CRISPR gene activation of mouse trans-2,3-enoyl-CoA reductase
Symbol:	Tocr
UniProt ID:	<u><a href="#">Q9CY27</a></u>
Synonyms:	2410016D23Rik; A230102P12Rik; AI173355; D17ErtD178e; Gpsn2; SC2
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

#### Product images:

