

Product datasheet for **GA211935**

Mouse Dlk2 activation kit by CRISPRa

Product data:

Product Type:	CRISPRa Kits
Description:	Dlk2 CRISPRa kit - CRISPR gene activation of mouse delta like non-canonical Notch ligand 2
Symbol:	Dlk2
Components:	GA211935G1, Dlk2 gRNA vector 1 in pCas-Guide-GFP-CRISPRa (GE100074) GA211935G2, Dlk2 gRNA vector 2 in pCas-Guide-GFP-CRISPRa (GE100074) GA211935G3, Dlk2 gRNA vector 3 in pCas-Guide-GFP-CRISPRa (GE100074) 1 CRISPRa-Enhancer vector, SKU GE100056 1 CRISPRa scramble vector, SKU GE100077
UniProt ID:	<u>Q8K1E3</u>
Synonyms:	AI413481; DLK-2; Egfl9
Format:	3 gRNAs (5ug each), 1 scramble ctrl (10ug) and 1 enhancer vector (10ug)

Product images:



[View online »](#)